

## General Immunofluorescence Staining Protocol using Directly Conjugated Antibodies

### 1. Test principle

Cells are stained using directly conjugated antibody for Immunofluorescence.

### 2. Specimen

Cells in suspension, from whole blood, bone marrow or cell culture

### 3. Materials and reagents

- 5ml Polypropylene test tubes (12x75 mm, round bottom)
- cooling centrifuge
- sterile-staining buffer (PBS, 2% FCS, 0.1 % azide) : The sodium azide assists in preventing capping and shedding or internalization of the antibody-antigen complex after the antibodies bind to the receptors.
- monoclonal antibody conjugated with fluorochrome
- Fc Blocker
- 1% paraformaldehyde in PBS
- flow cytometer

### 4. Controls

- Unstained cells
- Single stained controls in multicolor assays
- Isotypic IgG controls
- FMO control (= fluorescence minus one) in multicolor assays (Optional): Incubate with all colors except the one, you are interested in for that particular tube ⇒ negative control for color of interest

### 5. Procedure

1. Prepare single cell suspension and wash in staining buffer (PBS, 2% FCS, 0.1% azide).
2. Centrifuge (400 x g, 5 min, 4°C.), discard supernatant and resuspend by first flicking bottom of tube to break up pellet, and add staining buffer to  $1 \times 10^7$  cells/ml.
3. Aliquot 100  $\mu$ l of cells ( $10^6$  cells) into a 12 x 75 mm polypropylene tube.
4. **Blocking:** Add 5  $\mu$ l/tube of blocking antibody or serum (e.g. Fc Block). Pipetting up and down to mix and incubate for 2 min at room temp.
5. **Surface antibody staining:** Add monoclonal antibodies to the tube (turn off hood light for antibodies), pipetting up and down to mix and incubate for 30 min at 4°C **in the dark**.
6. Add 2 ml of staining buffer, mix well and centrifuge (400 x g, 5 min, 4°C).
7. Discard supernatant.
8. Wash again in 1 ml staining buffer and resuspend in 500  $\mu$ l staining buffer for FACS analysis.
9. Keep cells on ice prior to analysis.
10. Cells may be centrifuged and fixed in 1 ml of 1% paraformaldehyde (in PBS) at 4°C for analysis next day.

#### Note:

- Use buffers without Phenol Red
  - The blocking antibody step (4 and 5) is optional but should be included if cells express high levels of Fc receptors which will contribute to non-specific binding and background fluorescence.
  - You might need to adjust cell numbers, amount of antibody for your experiment.
-